

The Best Solution for Compression Therapy



“It’s a New Creation”

for the DVT prevention, edema therapy and quick recovery for affected area after surgery.

 **Doctor Supply**

Doctor Supply Quality Standards and Environment

Medical equipment manufacturing is a very competitive business that continuously demands advanced technology and functions.

Being a global supplier helps Doctor Supply to understand the complexity and competitive nature of the global medical equipment industry.

This understanding is the driving force behind our creative idea and innovative product development.

Also we believe that maintaining the high quality of our equipment through a highly standardized management system is the utmost priority.

Our various production systems are integrated as a single process and all of our business divisions are managed under the strict quality regulation of ISO 9001 and ISO 13485 which are international quality certifications. Besides, we have acquired CE Mark and FDA.

We have achieved the quality standardization for all our manufacturing facilities.



ISO 9001
ISO 13485



Product Lines

I Air Compression Line

- DOCTUS Air Cuff System
- Auto Com
- EZ-rap Revolution
- Thermo Cuff

III Compression Stocking Line

- CAP Stocking
- Health Support Stocking

II Elastic Compression Line

- EZ-rap



Doctor Supply Air Compression Line

Best solutions for DVT prevention and edema treatment after every surgery.

DOCTUS
Air Cuff System

DOCTUS
AUTOCom
Easy Control Automatic IPC System

EZrap
Revolution
Air Compression & PRICE Therapy

DOCTUS
Thermo Cuff
Best Solution for Edema Therapy after Surgery



DOCTUS

Air Cuff System



Best Solutions for DVT Prevention and Edema Treatment after Surgery.

DOCTUS Air Cuff System is a high-performance Intermittent Pneumatic Compression device and is capable of preventing DVT (Deep Vein Thrombosis), treating post-surgery edema.

- DVT Prevention
- Edema Therapy
- Quick Recovery for Surgery Effected Area

What is DVT (Deep Vein Thrombosis)?

DVT is the formation of a blood clot (thrombus) in a deep vein, most commonly of the calf and pelvis. Therefore, DVT and PE can arise mainly from knee or pelvis operation besides the neurological or other operations.

Healthy body suitably keeps the balance of the clotting system and the fibrinolytic system, but thrombus can arise when this balance breaks from several medical conditions that can lead to DVT, such as physical trauma, cancer, older age (the strongest risk factor).

Varicose vein, antiphospholipid syndrome, surgery, immobilization (as with bed-rest, orthopedic casts, or during long-haul flights), heavy smoking, obesity, certain drugs, and inborn tendencies form clots known as thrombophilia.

Anticoagulants

5.6%

Chance of Bleeding Complications

Intravenous Injections

95%

Numerous Undesirable Side Effects

IPC

0%

Chance of Bleeding Complications

Non-Invasive Therapy

95%

DVT Prevention

Practically **NO** Side Effects

VS.

Best effective and safe way to DVT prevention

"It's a New Creation"

Feature of DOCTUS Air Cuff System

DOCTUS Air Cuff System is composed of a controller which can operate and control all the functions of the system and Air Cuffs fitted around each of your legs or feet according to the symptom. The main body has two modes of DVT prophylaxis, edema treatment.

Each mode is designed to operate easily through the automatic air cuff sensory equipment.

Also, each patient to be served with various treatments in accordance with his or her status.

Feature of Air Cuff

Original air cuffs of DOCTUS Air Cuff System are designed for exclusive use by DOCTUS.

Totally, nine types of air cuffs of exclusive usage are broken down to DVT prevention only and edema treatment only.

All the Air Cuffs are produced with materials harmless to human body under our company's thoroughly strict sanitary and quality management.

Environment-friendly medical material



Adoption of ergonomic design



Sanitary care for patients by separate usage by individuals



Air Cuffs for DVT Prevention



Foot Cuff



Foot & Calf Cuff



Calf Cuff



Thigh Cuff

Air Cuffs for Edema Therapy



Boot Cuff



Full Leg Cuff



Arm Cuff




Full Arm Cuff

DOCTUS AutoCom

Easy Control Automatic IPC System



The best effective and easy homecare device for DVT prevention and edema treatment after surgery with the optimum air compression massage.



Various Availability

Applicable to all of Doctor Supplies' original air sleeves when connected to AutoCom.



DOCTUS Air Cuffs




EZrap Revolution



Thermo Cuff

"It's a New Creation"



Activity

- With mini type design ensuring portability
- Great usability in any place
- Long-time usability by built-in high-capacity battery

Economical

- Affordable high-performance air massage device
- Excellent economy ensured through sustainable use of a hospital-style air cuff at home
- Affordable Air Cuff for DVT prevention

Easily

- Easy Operation
- Optimized 5-step automatic modes for end-user
- Compatible with various conveniences
- With one-touch connector

Air Sleeves by Application

DVT Prevention

- After every abdominal surgery
- After every cancer surgery
- After total replacement hip & knee surgery
- Air compression Massage



Thigh Air Cuff



Calf Air Cuff



Foot Air Cuff



Foot & Calf Air Cuff

Edema Therapy


- After every abdominal surgery
- After every cancer surgery
- After total replacement hip & knee surgery
- Air compression Massage



Hip Thermo Cuff



Knee Thermo Cuff



Arm Air Cuff



Long Arm Air Cuff

Treatment of acute and chronic muscle

- After every abdominal surgery
- After every cancer surgery
- After total replacement hip & knee surgery
- Air compression Massage



Back Cuff



Shoulder Cuff



Knee / Elbow Cuff

EZrap Revolution

Air Compression & PRICE Therapy



The best solution for every muscular disease through air compression with cold and hot therapy system.

EZ-rap Revolution is the most efficient Cold & Hot Therapy product that makes use of air compression. Equipped with Air Cuff with Cold Pack used at the acute status (within 1 or 2 days after injury) and chronic status (in 2 days after injury), the product is designed to have maximum therapeutic effect by applying efficient air pressure on elbow, wrist, ankle, waist, shoulder, and higher sites while making use of air compression.

Particularly, the product is faithful to the orthopedic treatment principle of P.R.I.C.E (protection, rest, ice, compression, and elevation) and is easily put on and taken off.

And the product is available with hot or cold pack, which has shown the optimal effect in P.R.I.C.E on each joint of patients when moving.

Also, the product cools the heat produced within the joint of patient who moves, and the gel inside the Cold Pack is designed to stay immune to the gravity and spread evenly around the joint area, narrowing the gap between Air Cuff and the body part.

The most effective Cold/Hot Therapy is available by the ideal air compression that fits the actual status of the area.

"It's a New Creation"


Special Advantages and Features

- A line-up applicable to each part of the body.
- Designed to reflect ergonomic consideration, with the ideal shape for applying treatment using air pressure.
- Comfortable to put on; easy wearing adopted to enable a patient to put the product on and take it off by oneself.
- Use of harm-free medical materials provides protection from infection from outside.
- Cold Pack minimizes the gap between the Air Cuff and body part, enables the most effective cold or hot therapy with the fittest air pressure applied to sites of injury.



After Surgery

- It soothes tight joints and muscles and facilitates a person's movement with cold packs.



After Physical Therapy

- It reduces painful inflammation and bruises with cold packs.

Treating Pain, Stiffness, and Soreness

- It reduces swelling fast with compression and support.

What is P.R.I.C.E?

The P.R.I.C.E refers to the basic principle in orthopedic treatment. Our EZrap Revolution is designed in strict compliance with the orthopedic treatment principle.

- P** Protecting an affected area
- R** Providing rest to an affected area
- I** Applying an ice pack on an affected area
- C** Compression on an affected area
- E** Keep an affected area above the heart (Elevation)

Cold / Hot Pack (Soft Poly -Gel)



EZrap Revolution


Air Compression & PRICE Therapy

The EZrap Revolution can maximize the effect to when used connecting a AutoCom



Application for each body part

Air massage with cold & hot therapy are available at the same time.



DOCTUS
Thermo Cuff
Best Solution for Edema Therapy after Surgery



It is the most innovative and appropriate treatment for patients who have total knee replacement surgery and total hip replacement surgery.

Patients for total knee replacement surgery and total hip replacement surgery need a specialized therapy for edema therapy and DVT prevention after surgery.

Specifically, it can enhance the edema therapy and alleviate pain to a very appropriate level so that cold therapy and air massage can be applied at the same time.

Furthermore, Thermo Cuff has a unique design for edema therapy and DVT prevention for patients who under go total knee replacement surgery and total hip replacement surgery.

Also, Thermo Cuff guarantees a very outstanding effect with the optimum air massage through cold therapy for edema therapy after surgery for your home care.

Only One Solution!


for edema therapy and DVT prevention by optimum air massage and cold & hot therapy in a single package!



Air Massage



Cold Therapy



DOCTUS
ThermoCuff
 Best Solution for Edema Therapy after Surgery



Easy and consistent DVT prevention and edema treatment in your home following your discharge from hospital



Homecare

DVT prevention and all edema therapy available in your home.

Double Care

It increases therapeutic effect with optimum air massage and cold therapy when connected to AutoCom.



Knee Thermo Cuff



Hip Thermo Cuff

Feature and Indication

- Optimum ergonomic design for surgery patients.
- No skin irritation as it uses the finest medical fabric (Latex Free).
- Time and cost saved with convenient therapy and the ease in wearing the product.
- Capable of treating edema with support air compression.
- It quickly decreases edema and reduces pain.
- Applicable to breast cancer patients.
- Applicable to total knee and hip replacement surgery

Please use the product while it is connected to DOCTUS Air Cuff System or AutoCom.



Doctor Supply Elastic Compression Line

Best solution ensuring quick recovery for affected area with optimum elastic compression & cold therapy after surgery.

EZrap

Compression & Cold Therapy

- Abdominal Back Brace
- Ostomy Belt
- Shoulder Rap
- Chest Rap
- Elbow Rap
- Surgi Bra
- Hip Rap
- Bendable Wrist Splint
- Knee Rap / Knee Kool
- Knee 353



 Doctor Supply

EZrap

Compression & Cold Therapy



Best solution ensuring quick recovery for affected area with optimum elastic compression & cold therapy after surgery.

A variety of EZ-rap products which Doctor Supply produces provide exceptional effects in preventing discharge, blister after various kinds of surgery and reducing pain and edema. Particularly, products of EZ-rap prevent the rupture of a wound with intra-abdominal pressure and enable patients to resume walking earlier, having shown supreme effects in preventing post-operative lung-complications, intestinal adhesions, or abdominal hernia. (Especially, we commend that a patient after surgery puts on a stomachch band for six months.)

Blister prevention
On the wound

Surgical wound ruptured Prevention

Help in early ambulation

Reduce pain and edema

※ Weaknesses with existing cold therapies

- Frequent replacement is inevitable
- Ice therapy does not last long
- Considerable concern about blisters from serious condensation
- Ice causes pain to the skin
- Pressuring not available


"It's a New Creation"

Special Advantages and Features

- Prevents blisters around sites without taping for dressing purpose.
- Significant reduction of treatment time without extra elastic bandage.
- Provides parallel pain treatment with cold therapy function.
- Capable of treating edema with compression.
- Saves time and cost with convenient treatment through easy wearing.
- No irritation to skin as it uses the finest medical fabric (Latex Free)


Efficacy of Cold Therapy

Leads to contraction of muscle effective in reducing infection, edema, bleeding, and alleviation of pain.)




Tight compression achieved with high - quality elastic fabric.

Easy wearing by velcro tapping and no irritation to skin.



Guarantees mild cold effect with a soft multiple ice pack that does not freeze.

Smooth and constant pressure with special texture as it adopts magic velcro band that ensures fitness to patients' body.



EZrap

Compression & Cold Therapy

Abdominal Back Brace



◆ Feature and Benefit

- No need for taping with a tight pressure through flexible material.
- Absence of secondary injury or allergy by tight pressure around the dressing area.
- Lead to earlier recovery by enabling patients to work earlier without pain.
- Outstanding effect on pain reduction when deep breathing or coughing through a tight pressure.
- Superior effect on edema reduction by leading faster blood flow around the cut area.
- Cold Therapy function dramatically alleviates pain.

◆ Indication

- Back pain
- Stomach fat
- Abdominoplasty
- Liposuction surgery
- Pregnancy exercises
- Stomach pain

◆ Standard Specifications

Description	Length
Pediatric BACK RAP	21" (53 cm)
Small BACK RAP	21~27" (53~69 cm)
Medium BACK RAP	27~33" (69~84 cm)
Large BACK RAP	33~39" (84~99 cm)
x-Large BACK RAP	More than 39" (99 cm)

Abdominal Ostomy Belt



The specially designed and manufactured Ostomy Belt is for large bowel and hernia patients.

This particular belt is designed to provide a dedicated passage in urinary tract pocket with a view to removing discomfort in the tract pocket. And Ostomy belt is developed by a lot of clinicians based on the results of clinical demonstration.

◆ Indication

- Ostomy surgery



◆ Dedicated Passage for Urinary Tract Pocket after Ostomy Surgery.

"It's a New Creation"

Shoulder Rap Brace



◆ Feature and Benefit

- Similar effect with elastic band through adopting compression band easy to put on(reduction of edema by pressure), and fabulous blister preventing effect to the applied area.
- Cold Pack inside the SHOULDER RAP effective to reduction of strong pain and edema just after surgery guarantees close adhesion to the area.
- Lasting and strong cold therapy reduces infection, edema, and bleeding and alleviates pain dramatically.

◆ Indication

- After surgery of rotator cuff rupture
- After infection and pars Plana vitrectomy, other Shoulder abduction surgery

◆ Standard Specifications

Description	Length
Pediatric Shoulder Rap	27" (70cm)
Small Shoulder Rap	27~32" (70~80cm)
Medium Shoulder Rap	32~36" (80~90cm)
Large Shoulder Rap	36~43" (90~110cm)
x-Large Shoulder Rap	More than 43" (110cm)

Chest Rap



◆ Feature and Benefit

- Ergonomic design makes possible pressure for each parts of the body.
- Reusable after washing.
- Each tier of the product is less of distortion such as curling and wrinkling
- Minimum irritation to skin with the use of high-quality of cotton fabric.

◆ Indication

- After a variety of Open Heart surgery and lung surgery
- CABG (Coronary Artery Bypass Graft)
- After other chest surgery

◆ Standard Specifications

Description	Length
Pediatric Chest Rap	27" (70cm)
Small Chest Rap	27~32" (70~80cm)
Medium Chest Rap	32~36" (80~90cm)
Large Chest Rap	36~43" (90~110cm)
x-Large Chest Rap	More than 43" (110cm)

EZrap

Compression & Cold Therapy

Elbow Rap



◆ Feature and Benefit

- The product provides similar effects as with elastic band as it adopts a compression band easy to put on (reducing edema with pressure), and fabulous effect in preventing blisters in the area where it is applied.
- Cold Pack inside the ELBOW RAP, which is effective in reducing strong pain and edema immediately after surgery, guarantees close adhesion to the area.
- Lasting and strong cold therapy reduces infection, edema, and bleeding and alleviates pain dramatically.

◆ Indication

- Elbow arthroscopic surgery

◆ Standard Specifications

Description	Length
Pediatric Cool Rap	7" (17cm)
Small Cool Rap	7~9" (17~23cm)
Medium Cool Rap	9~11" (23~28cm)
Large Cool Rap	11~13" (28~33cm)
x-Large Cool Rap	More than 13" (33cm)

Surgi Bra



◆ Feature and Benefit

- Convenient to put on like usual garments by adopting ergonomic design
- Tightly fit to dressing area immovable after surgery
- Easy wearing by Velcro tapping and no irritation to skin

◆ Indication


- After surgery of Mammotomy
- After surgery of breast cancer

◆ Standard Specifications

Description	Length
Pediatric Surgi Bra	27" (70cm)
Small Surgi Bra	27~32" (70~80cm)
Medium Surgi Bra	32~36" (80~90cm)
Large Surgi Bra	36~43" (90~110cm)
x-Large Surgi Bra	More than 43" (110cm)

"It's a New Creation"

Hip Rap



◆ Feature and Benefit

- The product provides similar effects as with elastic band as it adopts a compression band easy to put on (reducing edema with pressure), and fabulous effect in preventing blisters in the area where it is applied.
- Cold Pack inside the HIP RAP, which is effective in reducing strong pain and edema immediately after surgery, guarantees close adhesion to the area .
- Lasting and strong cold therapy reduces infection, edema, and bleeding and alleviates pain dramatically.

◆ Indication

- THR (Total Hip Replacement Arthroplast)
- BPHA, Bipolar Hemiarthroplasty
- Fracture & Avascular necrosis Surgery
- Others Hip Surgery

◆ Standard Specifications

Description	Length
Pediatric Hip Rap	25" (64 cm)
Small Hip Rap	25~28" (64~71cm)
Medium Hip Rap	28~30" (71~76cm)
Large Hip Rap	30~35" (76~89 cm)
x-Large Hip Rap	More then 35" (89 cm)

Bendable Wrist Splint



◆ Feature and Benefit

- Tightly fix hand and wrist when inserting catheter and A Line.
- Keep needle from escaping or slipping out from infant or other unconscious patients.
- Capable of using as splint and very soft to wear.
- Various of size fit from infant to adult.
- Tailored application to patients by adopting elastic band.

◆ Indication

- Surgical patient of surgery and orthopedics (Including critical patient)

◆ Standard Specifications

Description	Length
Pediatric	for Infant
Small	for Children
Medium	for Adult

EZrap

Compression & Cold Therapy

Advantages for Knee Rap Line

- The product provides similar effects as with elastic band as it adopts a compression band easy to put on (reducing edema with pressure), and fabulous effect in preventing blisters in the area where it is applied.
- Cold Pack inside the KNEE RAP, which is effective in reducing strong pain and edema immediately after surgery, guarantees close adhesion to the area.
- Lasting and strong cold therapy reduces infection, edema, and bleeding and alleviates pain dramatically.
- Cold Therapy: leads to contraction of muscle effective to reduction of infection, edema, bleeding, and alleviation of pain.
- Ambulation is possible with cold pack.

◆ Indication

- After surgery of meniscal incision and repair
- After other surgery of knee joint arthroscopy
- After surgery of total knee joint arthroplasty
- After surgery of unicompartmental knee arthroplasty

◆ Standard Specifications

Description	Length
Pediatric Knee Rap	13" (33cm)
Small Knee Rap	13~15" (33~38cm)
Medium Knee Rap	15~17" (38~43cm)
Large Knee Rap	17~19" (43~48cm)
x-Large Knee Rap	More than 19" (48cm)



Knee Rap / Knee Kool

Prevention of dew condensation when applying EZ Pack pocket contained with cold therapy by using special fabric (two cold pack available)



Knee353

Three-way velcro type band adjustable degree of pressure and size to the fit of patients and symptoms. (three cold pack available)

Doctor Supply Compression Stocking Line

**The best high quality compression stocking
in the world-wide for the DVT prevention and
anti-pulmonary embolism .**

CAROLON
Health Support Stocking

Manufactured in CAROLON's state-of-the-art facility, stockings are tested in they special laboratory and certified for graduated compression by an industry leading external testing laboratory.

In utilizing the CAROLON system, each patient receives the highest quality protection against pulmonary embolism.

★ Doctor Supply is official agent of CAROLON for Asia market

CAP
Anti Embolism Stocking


HEALTH SUPPORT
Garadient Compression Stocking

 **Doctor Supply**



CAP

Anti Embolism Stocking



CAP anti embolism stocking

CAP stocking fit a larger portion of the patient population than any other anti-embolism stocking available today.

Knee length stockings are available in all standard size plus short, long and XL sizes.

This flexibility allows precision performance for both standard and difficult to fit patients.

During construction of the product, as many as fourteen different types of stitches and twelve different types of yarns are used.

This enables CAROLON to provide accurate controlled compression and patient comfort.

Feature and Indications

- COOL, LIGHTWEIGHT COMFORTABLE – Used to knit CAP stocking are strong yet lightweight for patient comfort.
- TWO WAY STRETCH FABRIC - Insures that the stocking remain in place, even during ambulation
- NON CONSTRICTIVE THIGH – Band helps hold stocking up
- FLOATING HEEL POCKET – Colored blue heel for easy recognition and accurate placement.
- CAP Stockings provide patients with certified graduated compression for the safe and effective prevention of DVT and Anti-pulmonary embolism.

Application of CAP



Knee High

Size	Ankle	Calf
S	17.5~20cm	22.5~33
M	20~22.5cm	30.5~38cm
L	22.5~25.5cm	35.5~43cm
XL	25.5~28cm	40.5~48cm



Thigh High

Size	Ankle	Calf	Thigh
S	17.5~20cm	22.5~33cm	40.5~51cm
M	20~22.5cm	30.5~38cm	48~61cm
L	22.5~25.5cm	35.5~43cm	58.5~68.5cm

Available

Small Short
Small Regular
Small Long

Medium Short
Medium Regular
Medium Long

Large Short
Large Regular
Large Long

X-Large Short
X-Large Regular
X-Large Long

XX-Large Short
XX-Large Regular
XX-Large Long

HEALTH SUPPORT

Gradient Compression Stocking

Comfortable and accurate of high functionality compression Stocking.

To ensure a high level of compliance, Health Support hosiery products utilize the most advanced yarn technologies and are constructed on state of the knitting equipment.


Products are developed with state of the art equipment to assure proper compression levels while providing maximum comfort and therapy.

Health Support hosiery products have been developed utilizing many of the advanced features used in the CAROLON anti-embolism stockings.

All Health Support hosiery products are machine washable and dryer safe.



Application of Health Support



Class 1
(20mmhg
~30mmhg)

1. Chronic leg fatigue
2. Mild venous insufficiency
3. Swelling of the ankles
4. Varicose veins

Available :
Knee, Sock, Full length thigh,
Panty hose, maternity panty hose,
Open / Closed toe

Class 2
(30mmhg
~40mmhg)

1. Chronic ankle swelling
2. Prophylaxis after vein stripping
3. Reversible lymph edema
4. Swelling after cast removal
5. Discomforts associated with phlebitis
6. Moderate to severe varicosities
7. Venous ulcerations

Available :
Knee, Sock, Full length thigh,
Panty hose, maternity panty hose,
Open / Closed Toe

Class 3
(40mmhg
~50mmhg)

1. Severe varicosities 2323
2. Severe edema
3. Severe venous ulceration
4. Lymph edema
5. Post thrombotic syndrome, chronic venous disease
6. Chronic venous insufficiency
7. Deep vein thrombosis

Available :
Knee, Full length thigh,

Doctor Supply Product Line-up

Air Compression Line

DOCTUS
Air Cuff System



Air Cuffs



Thigh Air Cuff



Calf Air Cuff



Foot Air Cuff



Foot & Calf
Air Cuff




Boot Air Cuff




Full Leg Air Cuff


DOCTUS
AUTOCom
Easy Control Automatic IPC System



Full Arm Cuff



Arm Air Cuff



Long Arm Air Cuff

"It's a New Creation"

EZrap Revolution Cuffs



Shoulder Cuff



Elbow Cuff



Back Cuff



Knee Cuff



Universal Cuff



Ankle Cuff



Hip Cuff

Thermo Cuffs



Knee Thermo Cuff



Hip Thermo Cuff

Elastic Compression with Cold Therapy Line

EZrap

Compression & Cold Therapy



Shoulder Rap



Abdominal Back Brace



Abdominal Ostomy Belt




Elbow Rap



Chest Rap



Surgi Bra



Hip Rap



Knee Rap
Knee Kool



Knee353

Compression Stocking Line

CAP

Anti Embolism Stocking



Knee High



Thigh High

HEALTH SUPPORT

Gradient Compression Stocking



Knee High



Thigh High



Panty Hose



Cushion Foot
Socks

EZpac

Optimum Cold & Hot Pack System




Cold Pack Solution

Guarantees mild cold effect with a soft multiple ice pack that does not freeze.

- No need ice frequent replacement.
- No pain to the skin by the ice
- Can be used for a long time
- Dew Condensation

◆ Application

- EZrap
- Thermo Cuff
- EZrap Revolution



Cold & Hot Pack Solution

EZpac is composed special poly gel and it is possible to choose cold pack or hot pack, depending on acute diseases and chronic diseases.

◆ Application

- EZrap
- Thermo Cuff
- EZrap Revolution





Doctor Supply Ltd.


39, Yangjae-daero 81gil, Gangdong-gu,
Seoul 134-847, Korea

Tel: +82-2-483-5834

Fax: +82-2-483-5836

Contact e-mail : pky@doctorsupply.co.kr

<http://www.doctorsupply.co.kr>



 **Doctor Supply**

"It's a New Creation"